

# SLC22A9 Antibody, HRP conjugated



PACO37699

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Sodium-independent organic anion transporter which exhibits high specificity for sulfated conjugates of xenobiotics and steroid hormones. It is also specifically activated by 3 to 5 carbons-containing short-chain fatty acid, /SCFAs, including propionate, butyrate and valerate. May operate the exchange of sulfated organic components against short-chain fatty acid, /SCFAs at the sinusoidal membrane of hepatocytes.

**Gene ID:**

SLC22A9

**Uniprot**

Q8IVM8

**Synonyms:**

Solute carrier family 22 member 9 (Organic anion transporter 7), SLC22A9, hOAT4 OAT7 UST3

**Immunogen:**

Recombinant Human Solute carrier family 22 member 9 protein (45-143AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A