

Abhd11 Antibody



PACO37722

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

mitochondrion, carboxylic ester hydrolase activity.

Gene ID:

Abhd11

Uniprot

Q8K4F5

Synonyms:

Protein ABHD11 (EC 3. -. -. -) (Alpha/beta hydrolase domain-containing protein 11) (Abhydrolase domain-containing protein 11) (Williams-Beuren syndrome chromosomal region 21 protein homolog), Abhd11, Wbscr21

Immunogen:

Recombinant Mouse Protein ABHD11 protein (1-307AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A