

# HERC6 Antibody, FITC conjugated



PACO37736

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates.

**Gene ID:**

HERC6

**Uniprot**

Q8IVU3

**Synonyms:**

Probable E3 ubiquitin-protein ligase HERC6 (EC 2.3.2.26) (HECT domain and RCC1-like domain-containing protein 6) (HECT-type E3 ubiquitin transferase HERC6), HERC6

**Immunogen:**

Recombinant Human Probable E3 ubiquitin-protein ligase HERC6 protein (764-1022AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A