

HERC6 Antibody, Biotin conjugated



PACO37737

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates.

Gene ID:

HERC6

Uniprot

Q8IVU3

Synonyms:

Probable E3 ubiquitin-protein ligase HERC6 (EC 2.3.2.26) (HECT domain and RCC1-like domain-containing protein 6) (HECT-type E3 ubiquitin transferase HERC6), HERC6

Immunogen:

Recombinant Human Probable E3 ubiquitin-protein ligase HERC6 protein (764-1022AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A