

Hmgcs1 Antibody, Biotin conjugated



PACO37765

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This enzyme condenses acetyl-CoA with acetoacetyl-CoA to form HMG-CoA, which is the substrate for HMG-CoA reductase.

Gene ID:

Hmgcs1

Uniprot

Q8JZK9

Synonyms:

Hydroxymethylglutaryl-CoA synthase, cytoplasmic (HMG-CoA synthase) (EC 2.3.3.10) (3-hydroxy-3-methylglutaryl coenzyme A synthase), Hmgcs1

Immunogen:

Recombinant Mouse Hydroxymethylglutaryl-CoA synthase, cytoplasmic protein (1-520AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A