

Dlk2 Antibody



PACO37770

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Regulates adipogenesis.

Gene ID:

Dlk2

Uniprot

Q8K1E3

Synonyms:

Protein delta homolog 2 (DLK-2) (Endothelial cell-specific protein S-1) (Epidermal growth factor-like protein 9) (EGF-like protein 9), Dlk2, Egf9


Immunogen:

Recombinant Mouse Protein delta homolog 2 protein (27-305AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: Dlk2 antibody at 5 μ g/ml + Mouse kidney tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 41, 46, 45 kDa. Observed band size: 41 kDa.