

KRTAP8-1 Antibody, Biotin conjugated



PACO37777

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

In the hair cortex, hair keratin intermediate filaments are embedded in an interfilamentous matrix, consisting of hair keratin-associated proteins (KRTAP), which are essential for the formation of a rigid and resistant hair shaft through their extensive disulfide bond cross-linking with abundant cysteine residues of hair keratins. The matrix proteins include the high-sulfur and high-glycine-tyrosine keratins.

Gene ID:

KRTAP8-1

Uniprot

Q8IUUC2

Synonyms:

Keratin-associated protein 8-1 (High glycine-tyrosine keratin-associated protein 8.1), KRTAP8-1, KAP8.1 KRTAP8.1

Immunogen:

Recombinant Human Keratin-associated protein 8-1 protein (1-63AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A