

TBC1D3 Antibody



PACO37802

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Acts as a GTPase activating protein for RAB5. Does not act on RAB4 or RAB11.

Gene ID:

TBC1D3

Uniprot

Q8IZP1

Synonyms:

TBC1 domain family member 3 (Prostate cancer gene 17 protein) (Protein TRE17-alpha) (Rab GTPase-activating protein PRC17), TBC1D3, PRC17 TBC1D3A

Immunogen:

Recombinant Human TBC1 domain family member 3 protein (1-291AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A