

# ANKMY2 Antibody



PACO37810

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May be involved in the trafficking of signaling proteins to the cilia.

**Gene ID:**

ANKMY2

**Uniprot**

Q8IV38

**Synonyms:**

Ankyrin repeat and MYND domain-containing protein 2, ANKMY2

**Immunogen:**


Recombinant Human Ankyrin repeat and MYND domain-containing protein 2 protein (1-280AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human liver cancer using PACO37810 at dilution of 1:100.