

IQCF2 Antibody



PACO37818

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

IQCF2

Uniprot

Q8IXL9

Synonyms:

IQ domain-containing protein F2, IQCF2

Immunogen:

Recombinant Human IQ domain-containing protein F2 protein (1-118AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A