

ARMC8 Antibody



PACO37826

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

ARMC8

Uniprot

Q8IUR7

Synonyms:

Armadillo repeat-containing protein 8, ARMC8

Immunogen:

Recombinant Human Armadillo repeat-containing protein 8 protein (1-240AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A