

PACO37842

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

CCDC36

Uniprot

Q8IYA8

Synonyms:

Interactor of HORMAD1 protein 1 (Cancer/testis antigen 74) (CT74) (Coiled-coil domain-containing protein 36), CCDC36, IHO1


Immunogen:

Recombinant Human Interactor of HORMAD1 protein 1 protein (11-592AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lung cancer using PACO37842 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO37842 at dilution of 1:100.