

# KLHDC8A Antibody, FITC conjugated



PACO37848

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Cul3-RING ubiquitin ligase complex, protein ubiquitination involved in ubiquitin-dependent protein catabolic process.

**Gene ID:**

KLHDC8A

**Uniprot**

Q8IYD2

**Synonyms:**

Kelch domain-containing protein 8A (Substitute for delta-EGFR expression 1) (S-delta-E1), KLHDC8A

**Immunogen:**

Recombinant Human Kelch domain-containing protein 8A protein (1-350AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A