

KLHDC8A Antibody, Biotin conjugated



PACO37849

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Cul3-RING ubiquitin ligase complex, protein ubiquitination involved in ubiquitin-dependent protein catabolic process.

Gene ID:

KLHDC8A

Uniprot

Q8IYD2

Synonyms:

Kelch domain-containing protein 8A (Substitute for delta-EGFR expression 1) (S-delta-E1), KLHDC8A

Immunogen:

Recombinant Human Kelch domain-containing protein 8A protein (1-350AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A