

KIAA0319L Antibody



PACO37850

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Possible role in axon guidance through interaction with RTN4R.

Gene ID:

KIAA0319L

Uniprot

Q8IZA0

Synonyms:

Dyslexia-associated protein KIAA0319-like protein (Adeno-associated virus receptor) (AAVR), KIAA0319L, AAVR KIAA1837


Immunogen:

Recombinant Human Dyslexia-associated protein KIAA0319-like protein (479-718AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO37850 at dilution of 1:100.



Immunofluorescent analysis of A549 cells using PACO37850 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).