

NOSTRIN Antibody, Biotin conjugated



PACO37857

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Multivalent adapter protein which may decrease NOS3 activity by inducing its translocation away from the plasma membrane.

Gene ID:

NOSTRIN

Uniprot

Q8IIV9

Synonyms:

Nostrin (BM247 homolog) (Nitric oxide synthase traffic inducer) (Nitric oxide synthase trafficker) (eNOS-trafficking inducer), NOSTRIN

Immunogen:

Recombinant Human Nostrin protein (1-322AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A