

# NKAIN4 Antibody, FITC conjugated



PACO37860

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

regulation of sodium ion transport.

**Gene ID:**

NKAIN4

**Uniprot**

Q8IVV8

**Synonyms:**

Sodium/potassium-transporting ATPase subunit beta-1-interacting protein 4 (Na(+)/K(+)-transporting ATPase subunit beta-1-interacting protein 4) (Protein FAM77A), NKAIN4, C20orf58 FAM77A

**Immunogen:**

Recombinant Human Sodium/potassium-transporting ATPase subunit beta-1-interacting protein 4 protein (83-150AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A