

ZC3H7A Antibody, FITC conjugated



PACO37868

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be a specific regulator of miRNA biogenesis. Binds to microRNAs MIR7-1, MIR16-2 and MIR29A hairpins recognizing the 3'-ATA(A/T)-5' motif in the apical loop.

Gene ID:

ZC3H7A

Uniprot

Q8IWR0

Synonyms:

Zinc finger CCCH domain-containing protein 7A, ZC3H7A, ZC3H7 ZC3HDC7

Immunogen:

Recombinant Human Zinc finger CCCH domain-containing protein 7A protein (1-167AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A