

CDH26 Antibody, HRP conjugated



PACO37879

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types.

Gene ID:

CDH26

Uniprot

Q8IXH8

Synonyms:

Cadherin-like protein 26 (Cadherin-like protein VR20), CDH26

Immunogen:

Recombinant Human Cadherin-like protein 26 protein (1-124AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A