

# CDH26 Antibody, FITC conjugated



PACO37880

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types.

**Gene ID:**

CDH26

**Uniprot**

Q8IXH8

**Synonyms:**

Cadherin-like protein 26 (Cadherin-like protein VR20), CDH26

**Immunogen:**

Recombinant Human Cadherin-like protein 26 protein (1-124AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A