

# COLGALT2 Antibody, HRP conjugated



PACO37891

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Has a beta-galactosyltransferase activity; transfers beta-galactose to hydroxylysine residues on collagen.

**Gene ID:**

COLGALT2

**Uniprot**

Q8IYK4

**Synonyms:**

Procollagen galactosyltransferase 2 (EC 2.4.1.50) (Collagen beta(1-O)galactosyltransferase 2) (ColGalT 2) (Glycosyltransferase 25 family member 2) (Hydroxylysine galactosyltransferase 2), COLGALT2, C1orf17 GLT25D2 KIAA0584

**Immunogen:**

Recombinant Human Procollagen galactosyltransferase 2 protein (351-626AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A