

COLGALT2 Antibody, Biotin conjugated



PACO37893

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Has a beta-galactosyltransferase activity; transfers beta-galactose to hydroxylysine residues on collagen.

Gene ID:

COLGALT2

Uniprot

Q8IYK4

Synonyms:

Procollagen galactosyltransferase 2 (EC 2.4.1.50) (Collagen beta(1-O)galactosyltransferase 2) (ColGalT 2) (Glycosyltransferase 25 family member 2) (Hydroxylysine galactosyltransferase 2), COLGALT2, C1orf17 GLT25D2 KIAA0584

Immunogen:

Recombinant Human Procollagen galactosyltransferase 2 protein (351-626AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A