

TCEAL8 Antibody



PACO37894

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

TCEAL8

Uniprot

Q8IYN2

Synonyms:

Transcription elongation factor A protein-like 8 (TCEA-like protein 8) (Transcription elongation factor S-II protein-like 8), TCEAL8


Immunogen:

Recombinant Human Transcription elongation factor A protein-like 8 protein (1-117AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human glioma using PACO37894 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human appendix tissue using PACO37894 at dilution of 1:100.