

TCEAL8 Antibody, Biotin conjugated



PACO37897

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in transcriptional regulation.

Gene ID:

TCEAL8

Uniprot

Q8IYN2

Synonyms:

Transcription elongation factor A protein-like 8 (TCEA-like protein 8) (Transcription elongation factor S-II protein-like 8), TCEAL8

Immunogen:

Recombinant Human Transcription elongation factor A protein-like 8 protein (1-117AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A