

CCDC155 Antibody



PACO37942

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in meiotic chromosome dynamics and homolog pairing.

Gene ID:

CCDC155

Uniprot

Q8N6L0

Synonyms:

Protein KASH5 (Coiled-coil domain-containing protein 155) (KASH domain-containing protein 5), CCDC155, KASH5

Immunogen:

Recombinant Human Protein KASH5 protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A