

LSMEM1 Antibody



PACO37954

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

cytoplasm.

Gene ID:

LSMEM1

Uniprot

Q8N8F7

Synonyms:

Leucine-rich single-pass membrane protein 1, LSMEM1, C7orf53


Immunogen:

Recombinant Human Leucine-rich single-pass membrane protein 1 protein (1-131AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human melanoma using PACO37954 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO37954 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).