

CLRN3 Antibody



PACO37974

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

extracellular exosome.

Gene ID:

CLRN3

Uniprot

Q8NCR9

Synonyms:

Clarín-3 (Transmembrane protein 12) (Usher syndrome type-3A-like protein 1), CLRN3, TMEM12 USH3AL1


Immunogen:

Recombinant Human Clarín-3 protein (39-91AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO37974 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).