

C3orf20 Antibody, HRP conjugated



PACO37979

Product Information

Size:	Protein Background:
50ug	cytoplasm.
Reactivity:	Gene ID:
Human	C3orf20
Source:	Uniprot
Rabbit	Q8ND61
Isotype:	Synonyms:
IgG	Uncharacterized protein C3orf20, C3orf20
Applications:	Immunogen:
ELISA	Recombinant Human Uncharacterized protein C3orf20 protein (239-493AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A