

OR10W1 Antibody



PACO37990

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Odorant receptor.

Gene ID:

OR10W1

Uniprot

Q8NGF6

Synonyms:

Olfactory receptor 10W1 (Olfactory receptor OR11-236), OR10W1, OR10W1P

Immunogen:

Recombinant Human Olfactory receptor 10W1 protein (152-188AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A