

C9orf43 Antibody



PACO38002

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

C9orf43

Uniprot

Q8TAL5

Synonyms:

Uncharacterized protein C9orf43, C9orf43


Immunogen:

Recombinant Human Uncharacterized protein C9orf43 protein (267-461AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human cervical cancer using PACO38002 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO38002 at dilution of 1:100.