

# FANK1 Antibody



PACO38018

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200,  
IF:1:50-1:500

**Protein Background:**

Plays an important in membrane trafficking through the secretory apparatus. Not involved in endocytic trafficking to lysosomes.

**Gene ID:**

FANK1

**Uniprot**

Q8TC84

**Synonyms:**

Fibronectin type 3 and ankyrin repeat domains protein 1, FANK1

**Immunogen:**


Recombinant Human Fibronectin type 3 and ankyrin repeat domains protein 1 protein (1-345AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human testis tissue using PACO38018 at dilution of 1:100.



Immunofluorescence staining of MCF-7 cells with PACO38018 at 1:333, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).