

MRPL30 Antibody, Biotin conjugated



PACO38025

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

mitochondrial inner membrane, mitochondrion, mitochondrial translation, mitochondrial translational elongation, mitochondrial translational initiation, mitochondrial translational termination, organelle organization.

Gene ID:

MRPL30

Uniprot

Q8TCC3

Synonyms:

39S ribosomal protein L30, mitochondrial (L30mt) (MRP-L30) (39S ribosomal protein L28, mitochondrial) (L28mt) (MRP-L28) (Mitochondrial large ribosomal subunit protein uL30m), MRPL30, MRPL28 RPML28

Immunogen:

Recombinant Human 39S ribosomal protein L30, mitochondrial protein (35-161AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A