

SFXN5 Antibody



PACO38026

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transports citrate. Potential iron transporter.

Gene ID:

SFXN5

Uniprot

Q8TD22

Synonyms:

Sideroflexin-5, SFXN5

Immunogen:

Recombinant Human Sideroflexin-5 protein (1-102AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A