

# AEN Antibody, HRP conjugated



PACO38047

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Exonuclease with activity against single- and double-stranded DNA and RNA. Mediates p53-induced apoptosis. When induced by p53 following DNA damage, digests double-stranded DNA to form single-stranded DNA and amplifies DNA damage signals, leading to enhancement of apoptosis.

**Gene ID:**

AEN

**Uniprot**

Q8WTP8

**Synonyms:**

Apoptosis-enhancing nuclease (EC 3.1. -. -) (Interferon-stimulated 20 kDa exonuclease-like 1), AEN, ISG20L1

**Immunogen:**

Recombinant Human Apoptosis-enhancing nuclease protein (117-325AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A