

AEN Antibody, Biotin conjugated



PACO38049

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Exonuclease with activity against single- and double-stranded DNA and RNA. Mediates p53-induced apoptosis. When induced by p53 following DNA damage, digests double-stranded DNA to form single-stranded DNA and amplifies DNA damage signals, leading to enhancement of apoptosis.

Gene ID:

AEN

Uniprot

Q8WTP8

Synonyms:

Apoptosis-enhancing nuclease (EC 3.1. -. -) (Interferon-stimulated 20 kDa exonuclease-like 1), AEN, ISG20L1

Immunogen:

Recombinant Human Apoptosis-enhancing nuclease protein (117-325AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A