

TMEM255B Antibody



PACO38058

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Gene ID:

TMEM255B

Uniprot

Q8WV15

Synonyms:

Transmembrane protein 255B (Protein FAM70B), TMEM255B, FAM70B


Immunogen:

Recombinant Human Transmembrane protein 255B protein (221-356AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunofluorescent analysis of MCF-7 cells using PACO38058 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO38058 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO38058 at dilution of 1:100.