

EEF1AKMT1 Antibody, Biotin conjugated



PACO38073

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

S-adenosyl-L-methionine-dependent protein-lysine N-methyltransferase that methylates elongation factor 1-alpha.

Gene ID:

EEF1AKMT1

Uniprot

Q8WVE0

Synonyms:

EEF1A lysine methyltransferase 1 (EC 2.1.1. -) (N(6)-adenine-specific DNA methyltransferase 2) (Protein-lysine N-methyltransferase N6AMT2), EEF1AKMT1, N6AMT2

Immunogen:

Recombinant Human EEF1A lysine methyltransferase 1 protein (2-214AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A