

DYDC1 Antibody



PACO38078

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Plays a crucial role during acrosome biogenesis.

Gene ID:

DYDC1

Uniprot

Q8WWB3

Synonyms:

DPY30 domain-containing protein 1, DYDC1, DPY30D1


Immunogen:

Recombinant Human DPY30 domain-containing protein 1 protein (1-177AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: DYDC1 antibody at 4 μ g/ml + HeLa whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 21, 22 kDa. Observed band size: 21 kDa.