

LPGAT1 Antibody, FITC conjugated



PACO38092

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Lysophosphatidylglycerol (LPG) specific acyltransferase that recognizes various acyl-CoAs and LPGs as substrates but demonstrates a clear preference for long chain saturated fatty acyl-CoAs and oleoyl-CoA as acyl donors. Prefers oleoyl-LPG over palmitoyl-LPG as an acyl receptor and oleoyl-CoA over lauroyl-CoA as an acyl donor.

Gene ID:

LPGAT1

Uniprot

Q92604

Synonyms:

Acyl-CoA: lysophosphatidylglycerol acyltransferase 1 (EC 2.3.1. -), LPGAT1, FAM34A KIAA0205

Immunogen:

Recombinant Human Acyl-CoA: lysophosphatidylglycerol acyltransferase 1 protein (151-341AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A