

BORCS5 Antibody



PACO38106

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

BORCS5

Uniprot

Q969J3

Synonyms:

BLOC-1-related complex subunit 5 (Loss of heterozygosity 12 chromosomal region 1) (Myristoylated lysosomal protein) (Myrlysin), BORCS5, LOH12CR1


Immunogen:

Recombinant Human BLOC-1-related complex subunit 5 protein (1-196AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO38106 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO38106 at dilution of 1:100.