

MAGOHB Antibody, Biotin conjugated



PACO38113

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in mRNA splicing and in the nonsense-mediated decay (NMD) pathway.

Gene ID:

MAGOHB

Uniprot

Q96A72

Synonyms:

Protein mago nashi homolog 2, MAGOHB, MAGOH2

Immunogen:

Recombinant Human Protein mago nashi homolog 2 protein (2-148AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A