

S100PBP Antibody



PACO38114

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

nucleus, calcium-dependent protein binding.

Gene ID:

S100PBP

Uniprot

Q96BU1

Synonyms:

S100P-binding protein (S100P-binding protein Riken), S100PBP, S100PBPR


Immunogen:

Recombinant Human S100P-binding protein (276-408AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lung cancer using PACO38114 at dilution of 1:100.