

RAB39B Antibody, HRP conjugated



PACO38131

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in vesicular trafficking. Plays a role in synapse formation. May regulate the homeostasis of SNCA/alpha-synuclein. Together with PICK1 proposed to ensure selectively GRIA2 exit from the endoplasmic reticulum to the Golgi and to regulate AMPAR composition at the post-synapses and thus synaptic transmission.

Gene ID:

RAB39B

Uniprot

Q96DA2

Synonyms:

Ras-related protein Rab-39B, RAB39B

Immunogen:

Recombinant Human Ras-related protein Rab-39B protein (1-213AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A