

# DPH3 Antibody, HRP conjugated



PACO38147

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Essential for the first step in the synthesis of diphthamide, a post-translational modification of histidine which occurs in elongation factor 2 (EEF2) and which can be ADP-ribosylated by diphtheria toxin and by Pseudomonas exotoxin A (Eta).

**Gene ID:**

DPH3

**Uniprot**

Q96FX2

**Synonyms:**

DPH3 homolog (CSL-type zinc finger-containing protein 2) (DelGEF-interacting protein 1) (DelGIP1), DPH3, DESR1 ZCSL2

**Immunogen:**

Recombinant Human DPH3 homolog protein (1-82AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A