

DPH3 Antibody, Biotin conjugated



PACO38149

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Essential for the first step in the synthesis of diphthamide, a post-translational modification of histidine which occurs in elongation factor 2 (EEF2) and which can be ADP-ribosylated by diphtheria toxin and by Pseudomonas exotoxin A (Eta).

Gene ID:

DPH3

Uniprot

Q96FX2

Synonyms:

DPH3 homolog (CSL-type zinc finger-containing protein 2) (DelGEF-interacting protein 1) (DelGIP1), DPH3, DESR1 ZCSL2

Immunogen:

Recombinant Human DPH3 homolog protein (1-82AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A