

HHIPL1 Antibody, Biotin conjugated



PACO38169

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	HHIPL1
Source:	Uniprot
Rabbit	Q96JK4
Isotype:	Synonyms:
IgG	HHIP-like protein 1, HHIPL1, HHIP2 KIAA1822
Applications:	Immunogen:
ELISA	Recombinant Human HHIP-like protein 1 protein (401-608AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A