

MS4A4A Antibody, FITC conjugated



PACO38172

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in signal transduction as a component of a multimeric receptor complex.

Gene ID:

MS4A4A

Uniprot

Q96JQ5

Synonyms:

Membrane-spanning 4-domains subfamily A member 4A (CD20 antigen-like 1) (Four-span transmembrane protein 1), MS4A4A, 4SPAN1 CD20L1 MS4A4

Immunogen:

Recombinant Human Membrane-spanning 4-domains subfamily A member 4A protein (1-64AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A