

CAPNS2 Antibody, HRP conjugated



PACO38179

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Calcium-regulated non-lysosomal thiol-protease which catalyzes limited proteolysis of substrates involved in cytoskeletal remodeling and signal transduction. This small subunit may act as a tissue-specific chaperone of the large subunit, possibly by helping it fold into its correct conformation for activity.

Gene ID:

CAPNS2

Uniprot

Q96L46

Synonyms:

Calpain small subunit 2 (CSS2) (Calcium-dependent protease small subunit 2), CAPNS2

Immunogen:

Recombinant Human Calpain small subunit 2 protein (1-248AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A