

# SENP8 Antibody, FITC conjugated



PACO38184

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Protease that catalyzes two essential functions in the NEDD8 pathway: processing of full-length NEDD8 to its mature form and deconjugation of NEDD8 from targeted proteins such as cullins or p53.

**Gene ID:**

SENP8

**Uniprot**

Q96LD8

**Synonyms:**

Sentrin-specific protease 8 (EC 3.4.22. -) (Deneddylase-1) (NEDD8-specific protease 1) (Protease, cysteine 2) (Sentrin/SUMO-specific protease SENP8), SENP8, DEN1 NEDP1 PRSC2

**Immunogen:**

Recombinant Human Sentrin-specific protease 8 protein (1-212AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A