

SERPINB12 Antibody, HRP conjugated



PACO38191

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Inhibits trypsin and plasmin, but not thrombin, coagulation factor Xa, or urokinase-type plasminogen activator.

Gene ID:

SERPINB12

Uniprot

Q96P63

Synonyms:

Serpin B12, SERPINB12

Immunogen:

Recombinant Human Serpin B12 protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A